

### Smart Clinical Board TKG-107-TOUCH-IP68-WHITE







### **PRODUCT FEATURES:**

- Keyboard with silicone enclosure
- Completely waterproof and dustproof
- Can be fully disinfected and cleaned
- Integrated touchpad

- Locking key to prevent any incorrect commands being triggered
- Flat key caps are easy to clean
- High quality manufacture
- 3 years warranty on all electrical parts

#### **TECHNICAL DATA**

Power supply:

Number of keys:	107
Switch technology:	carbon contacts on gold
Shift force:	1.5 N
Shift path:	1.3 mm
Service life:	each key can be used approx.
	1 million times
Protection class:	IP68
Enclosure material:	silicone
Enclosure colour:	white
Interface:	USB

100 mA

Cable length:	2 m
Weight:	ca. 913 g
Operating temperature:	0 °C to + 70 °C
Storage temperature:	-30 °C to + 70 °C
Dimensions:	377 x 157 x 15 mm
TECHNICAL DATA FOR TOUC Operating principle: Dimensions:	CHPAD capacitive 60 x 42 mm
Documentation:	CE declaration, IP test, guarantee certificate, EMC test


Cat. no.	Product description	Keyboard layout
KG22204	TKG-107-TOUCH-IP68-WHITE-USB-DE	DE
KG22205	TKG-107-TOUCH-IP68-WHITE-USB-US/EU	US/EU
Other layouts, configurations and interfaces are available on request		

# www.ipc2u.ru www.ipc2u.de www.ipc2u.com

## Date. 12/2015 Rev.01







## **1. KEY FUNCTIONS:**



#### 2. SPECIAL FUNCTIONS:

Clean mode:	To activate or deactivate the clean mode, hold this key ad down for 2 seconds. This deactivates the key function and the touchpad. You can then clean the keyboard without triggering any incorrect commands.
*Touchpad-tapping:	Touchpad tapping allows you to use different functions by typing the touchpad surface once or several times (e.g. in word processing: single tap to position the cursor, double tap to select paragraphs)

You can find operating instructions and hints on cleaning and disinfecting at:

www.gett-group-service.com

GETT Gerätetechnik GmbH Mittlerer Ring 1 08233 Treuen Germany Telefon/Phone: Telefax/Fax: E-Mail: +49 37468 660-0 +49 37468 660-66 info@gett.de

gett-group.com